BATHROOM APPLIANCES



Model		BAS16-A	BA16-C
	Rear cleansing	1 Minute Max, 0.6 LPM (0.16 GPM)	
Operating Time	Front cleansing	Self-cleaning nozzle	
	Warm Air Drying	2 Minutes	
Control	Move	0	
	Nozzle Position	5 Levels	
	Water pressure	3 Levels	
	Water Temperature	4 Levels (Off : Room temperature, L : 32°C, M : 35°C, H : 38°C)	
	Air Dryer Temperature	3 Levels (L : 32°C, M : 35°C, H : 38°C)	
	Heated Seat	4 Levels (Off : Room temperature, L : 31°C, M : 34°C, H : 37°C)	
Sterilization Bowl		0	Χ
Air Pump		0	Χ
Nozzle Cleaning		0	
Filters		MF Filter	
Power Consumption		1,050 W (Max)	
Operating Noise		37 ~ 48 dB	
Net Weight		5.3 kg (11.7 lb)	
Dimensions (W X D X H)		190 190-110	494 X 534 X 150 mm (19.4 X 21.0 X 5.9 inch)

About COWAY

Made with Coway's exceptional core technology, our highly competitive products are popular in more than 80 countries around the world. As the world's largest specialist in water filtration appliances and leading edge household well-being electronics, we have developed a product line that meets the diversified needs of global consumers through relentless research in environment-friendly products across a broad range of categories.

Increasingly recognized as a world leader in innovative household health / lifestyle technology, Coway's awards and special recognitions continue to accumulate.

AWARDS











































Advanced Sterilization Bidet - BAS16





BAS16-A/C STERILIZATION BOWL BIDET





ENJOY IMPROVED HEALTH, COMFORT, AND IMPECCABLE CLEANLINESS

These microbes breed in the toilet bowl and contaminate the air after flushing. To prevent this, the BAS-16 bidet sterilizes microbes to provide a hygienic and clean washroom environment.

ELECTROLYSIS STERILIZING CARE SYSTEM *For BAS16-A only


Chemical-free & Eco-Friendly technology sterilizes your toilet from microbes such as E. Coli, Faecalis and other microorganisms with oxidants that do not harm the environment.




The sterilizing water is generated by a

simple push of the bowl sterilize button. * WHAT IS COWAY STERILIZATION IONIZER (CSI)?

Coway introduces an ionized-water (Mixed Oxidants) sterilization technology to sterilize the toilet bowl without any chemicals.



The sterilizing water comes out from the nozzles and is sprayed into the bowl.



The entire toilet bowl is sterilized for 15 seconds.



COWAY STERILIZATION BIDET VS. GENERAL BIDET BAS16-A

Eliminates over 99.9% of microbes such as E. Coli, Faecalis.





SUPERIOR HYGIENE SYSTEMS Stainless Steel Twin nozzles and an anti-bacteria seat ensure daily cleanness.



Nozzle Door

The newly applied nozzle door protects the nozzle part that is the most frequently exposed to contaminants.



Detachable Nozzle Duct

removable nozzle duct and can easily be replaced and ed in a single touch.



Stainless steel twin nozzle


The hardened stainless steel prevents bacterial contamination and foreign matters from building up.



leaning nozzle

e automatically creates ensure its cleanliness g the bidet.





ELECTRONIC BIDET SYSTEMS Adaptable functions, ease of attachment



Air* water stream

An air pump injects bubbles into the water stream fo



Easy detachment

With just one touch a user can remove and replace



Heated seat

The temperature control device maintains comfort.



Universal Bracket



Dry control

After the cleansing is complete, the temperature adjustable dry air guarantees a refreshing finish.



Anti-bacteria Seat

Anti-bacteria seats restrict bacteria propagation to provide consistent